



Princeton Instruments

Spec-10:100B

Cryogenic Cooling
1340 x 100 imaging array
20 x 20- μ m pixels

The Princeton Instruments Spec-10:100B utilizes a high-performance, back-illuminated, spectroscopic-format CCD designed exclusively for Roper Scientific®. The 1340 x 100 imaging array is ideal for spectroscopy applications, providing small chip height for fast spectral rates along with 27-mm spectral coverage. This detector delivers much higher resolution and sensitivity than industry-standard "1024 pixel" sensors. A second exclusive feature is the integration of software-selectable amplifiers that offer an easy choice of high sensitivity or high signal-to-noise ratio (SNR). Liquid nitrogen cooling eliminates dark noise, even for long exposure times.

Features	Benefits
1340 x 100 imaging array	Exclusive feature that provides superior resolution over industry-standard "1024 pixel" format
Exclusive CCD architecture	Provides industry's lowest-noise CCD system
20 x 20- μ m pixels	Optimum pixel size for full well and high resolution
2-mm chip height	Ideal for rapid spectral acquisition
Software-selectable amplifiers	Exclusive feature provides choice of superior sensitivity or superior SNR
Back-illuminated CCD	High quantum efficiency
Cryogenic cooling	Eliminates noise attributable to dark current, even for long exposure times
Standard spectrometer interface	Will interface with most spectrometers
Dual-digitizer option	Multiple-speed digitization allows complete freedom to select between "slow operation" for low noise and highest SNR or "fast operation" for rapid image acquisition
"USB 2.0 interface" configuration	Seamless, plug-and-play connection to PC notebooks and desktops Easy OEM integration
"PCI interface" configuration	Industry standard for fast, reliable data transfer
WinSpec and PVCAM®	Offers easy-yet-sophisticated Windows® GUI controls Automates data acquisition, analysis, and display
Linux® drivers and SITK™ plug-in for National Instruments' LabVIEW™	Extends system utility

Specifications

CCD image sensor	Princeton Instruments exclusive; scientific grade 1; MPP; back illumination; available with Unichrome UV-enhancement coating				
CCD format	1340 x 100 imaging pixels; 20 x 20- μ m pixels; 100% fill factor; 26.8 x 2.0-mm imaging area				
	Minimum		Typical		Maximum
Dual digitizers: system read noise			2.8 e- rms		3 e- rms
@ 50-kHz digitization			3.5 e- rms		5 e- rms
@ 100-kHz digitization			5 e- rms		6 e- rms
@ 200-kHz digitization			6 e- rms		8 e- rms
@ 500-kHz digitization			8 e- rms		10 e- rms
@ 1-MHz digitization			13 e- rms		16 e- rms
Single digitizer: system read noise			2.8 e- rms		3 e- rms
@ 50-kHz digitization			3.5 e- rms		5 e- rms
@ 100-kHz digitization			8 e- rms		10 e- rms
@ 1-MHz digitization					
Spectrometric well capacity	high sensitivity	high capacity	high sensitivity	high capacity	
	250 ke-	800 ke-	300 ke-	1 Me-	
Deepest cooling temperature			-120°C		
Dark current @ -120°C operation			0.3 e-/p/hr		1 e-/p/hr
Software-selectable gains	1/2x, 1x, 2x (high-sensitivity mode) 1x, 2x, 4x (high-capacity mode)				
Dynamic range	16 bits				
Nonlinearity	<2%				
Vertical shift time	30 μ s				
Spectral rates*					
@ 100-kHz digitization	60 Hz				
@ 1-MHz digitization	230 Hz				
@ 2-MHz digitization	270 Hz				
Operating temperature	-70 to -120°C				
Thermostating precision	\pm 0.05°C across entire temperature range				

SPEC-10
PRINCETON INSTRUMENTS



Note: Specifications are subject to change.

*Spectral rates have been measured with 100 rows vertically binned.

When you're **SERIOUS** about high-performance spectroscopy...

www.roperscientific.com

USA East Coast 609.587.9797 **USA West Coast** 520.889.9933 **Benelux** +31.347.324989
France +33.160.86.03.65 **Germany** +49.89.660.779.3 **Japan** +81.43.274.8022

Spec-10:100B (LN) Rev D1

Princeton Instruments Spec-10:100B (LN)

