

HTS Lens Spectrograph



The Princeton Instruments / Acton Research HTS Lens Spectrograph is an f/1.8 imaging spectrograph optimized for fiberoptic light collection. The HTS provides superior imaging performance in a compact, rugged package that includes an internal shutter, manually controlled wavelength adjustment, and a manually adjustable bilateral entrance slit. Two interchangeable gratings are available, allowing operation over a range of approximately 400 nm to 850 nm. A set of wavelength calibration tables is provided with each system, facilitating simple integration into your experimental setup. This spectrograph is supported by the industry-standard SpectraSense® and WinSpec software packages. The HTS can be integrated easily with a wide range of accessories, several of which include a notch filter chamber or fiberoptic input adapter.

Applications: Raman, luminescence, multistripe emission spectroscopy

Features	Benefits
Superior spatial imaging	Well suited for multiple-input measurements
Fast f/1.8 optical spectrograph	Ideally matched for fiberoptic input
Compact, rugged design	Simple integration into applications where space is at a premium
Interchangeable grating turret	Allows you to optimize for the best dispersion
Less than 100-µm astigmatism	Image height remains constant across the tangential focal plane
Tunable center wavelengths	Usable range from 400 nm to 850 nm
Renowned WinSpec and SpectraSense® software	Offers easy, yet sophisticated Windows® GUI controls Automates data acquisition, analysis, and display
PICAM for VB, C, C++, and Scientific Toolkit for LabVIEW™	Respected application programming interface provides a universal interface for all Princeton Instruments / Acton Research hardware
Wide range of accessories available	Includes fiber adapters, sample chambers, and notch filters

HTS Lens Spectrograph Rev A0

1


HTS Lens Spectrograph Specifications

Focal length	85 mm
Aperture ratio	f/1.8
Nominal dispersion 600 g/mm 830 g/mm	10 nm/mm 6 nm/mm
Spectral resolution (10-µm slit) 600 g/mm 830 g/mm	1.0 nm 0.5 nm
Spatial performance (astigmatism)	<100 µm
Focal plane size	28 mm x 8 mm
Grating size	32 mm x 32 mm
Grating mount	Interchangeable grating holder
Wavelength adjustment	Micrometer controlled, manual

All specifications are subject to change without notice.

PI*MAX and Spec-10 are registered trademarks of Roper Scientific, Inc. PhotonMAX and PIXIS are trademarks of Roper Scientific, Inc. SpectraSense is a registered trademark of Acton Research Corporation. LabVIEW is a trademark of National Instruments Corporation. Windows is a registered trademark of Microsoft Corporation in the United States and other countries. Other brand and product names are the trademarks or registered trademarks of their respective owners and manufacturers.

Outline Drawings



200

400

600

800

Wavelength (nm)

Diffraction Gratings Range of operation Wavelength coverage Part number Groove density Blaze wavelength Dispersion (1-inch focal plane) HTS1-060-VIS 600 g/mm 500 nm 350 – 750 nm 10 nm/mm 254 nm HTS1-083-VIS 830 g/mm 820 nm 500 - 900 nm 6 nm/mm 152 nm **HTS Gratings** 100 80 Efficiency (%) 600-g/mm ruled grating, blazed at 500 nm 60 40 830-g/mm ruled grating, blazed at 820 nm 20

1000

1200

1400

				Accessorie	es	
		Part number		Description		
Raman notch filter chamber		NFC-446-040			Raman notch filter assembly, dual lens, slit mount (requires notch filter)	
Туре	Part nun	nber	Wavelength	Laser attenuation (optical density)	Spectral bandwidth @ 50% transmission	
Notch	NFH-N-514.5 NFH-N-532 NFH-N-632.8 NFH-N-785		514.5 nm 532 nm 632.8 nm 785 nm	>4.0	700 nm (20 cm ⁻¹)	
Notch Plus	NFH-NP-514.5 NFH-NP-532 NFH-NP-632.8 NFH-NP-785		514.5 nm 532 nm 632.8 nm 785 nm	>6.0	700 nm (20 cm ⁻¹)	
Super Notch	NFH-SN-S NFH-SN-S	532	514.5 nm 532 nm 632.8 nm	>4.0	350 nm (10 cm ⁻¹)	

785 nm

514.5 nm

632.8 nm

532 nm

785 nm



NFH-SN-785

NFH-SNP-514.5

NFH-SNP-632.8

NFH-SNP-532

NFH-SNP-785

Super Notch Plus

www.piacton.com

>6.0

email: moreinfo@piacton.com USA +1.609.587.9797 | Benelux +31 (347) 324989 France +33 (1) 60.86.03.65 | Germany +49 (0) 89.660.779.3 Japan +81.3.5639.2741 | UK +44 (0) 28.383101.71

350 nm (10 cm⁻¹)